COMBINES FOUR DIFFERENT TECHNOLOGIES RF+Vacuum+IR+ **Ultrasound Cavitation 4 in 1** Skin Tightening Fat Reduction Cellulitis Treatment Body Contouring Face Lifting Wrinkles Removal

12.1inch high-resolution capacitive screen (touchable)





JAPAN SMC SOLENOID VALVE

High speed, long life, no noise



ITALY IMPORTED AIR PUMP


Large flow, low noise



AMERICA COHERENT IR LASER DIODE

High power, long life





Intelligent Operating System

• The handle can be automatically recognized, injection molding process, lightweight and convenient


Japan Imported Triple Frequency RF Power Supply



• It adopts triple frequency RF power supply imported from Japan, which can output three different frequencies of 1M, 5M and 8M respectively, and three different penetration depths. Customers can choose different frequencies for treatment according to different treatment items.










- 1M RF for body slimming, weight loss and body shaping
- 5M RF for deep wrinkles wrinkle removal, skin tightening and lifting
- 8M RF for superficial wrinkle removal and skin tightening

Vacuum RF Handle







Combination of Multi-polar RF, Vacuum & IR makes it possible to deliver RF energy to the desired target area effectively and precisely. Both face and body treatment is possible, offering a wide range of benefits such as anti-wrinkle and skin regeneration effect.




Milti-Polar RF Handle






The combination of Multi-polar RF and IR enables accurate and effective RF energy transfer while improving skin elasticity. 3 size of handpieces allow various area treatment from face to whole body.



Ultrasound Cavation Handle





- 60Khz ultrasonic wave, directly reach the fat layer, the fat reduction effect is remarkable.
- 80mm diameter ultra-large ultrasonic array, provides better whole-body fat reduction effect, and faster treatment speed.

Before And After







After



Before



After





Before



After



Before



After



Before



After



Before



After



Before



After



Before



After

Specifications

Voltage: 220V/110V

Power: 500W

Vacuum: 82KP

RF frequency: 1MHz/5MHz/8MHz

Output mode: Continuous output

RF power: 300W

Screen: 12.1 inch Touch screen

Handle quantity: 7 pieces

