

SMAS FOCUS V

/ Skin Tightening · Face Lifting · Body Slimming · Vaginal Tightening /


5D HIFU



One Machine, Multi-Functional

Only need to replace the handle, you can upgrade the product function at any time.


How it works



HIFU-powered transducers deliver concentrated waves of ultrasound into multiple layers of skin and body tissues without making direct contact with surrounding areas.

Transferring stabilized temperatures of 65~75°C into the targeted zones, the SMASFOCUS V aims to induce a process of neocollagenesis in the face while tightening body tissues for contouring purposes

FOR LIFTING & TIGHTENING



1 The dermis layer of the skin naturally synthesizes less collagen over the course of time.



3 Stimulated by the focused energy, skin collagen starts to regeneration.




2 Focused ultrasound waves deliver energy directly within the skin layer causing a shrinking effect.



4 Through the regeneration of collagen, wrinkles are eliminated in the process of both firming and tightening the skin.


FOR SLIMMING & CONTOURING



1 The dermis layer of the skin naturally synthesizes less collagen over the course of time.



3 As a natural healing response of the body the cells regenerate joining together to fill the gap.




2 The SMASFOCUS V energy destroys the fat cells that have been targeted.



4 In result, the layer of fat naturally exits the body delivering a smooth, contoured body line.

FOR VAGINAL TIGHTENING





Multiple choices of handles



Vmax single dot handle



1.5mm, 3.0mm, 4.5mm



2D hifu handle



1.5mm, 3.0mm, 4.5mm



3D hifu handle



1.5mm, 3.0mm, 4.5mm



6.0mm, 8.0mm, 10.0mm, 13.0mm, 16.0mm



Liposonix handle



8.0mm, 13mm



Vaginal handle



3.0mm, 4.5mm

Eyebrow, Wrinkles around eyes,
Wrinkles around month


Check, Forehead, Double chin,
Neck wrinkles, V-line forming

Check, Forehead, Double chin,
Neck wrinkles, V-line forming

Abdomen, Inner thighs, Outer
thighs, Hips & Buttocks, Inner
arms, Breast


Body slimming and fat reduction

Vaginal tightening




Easy to upgrade


Each machine is pre-loaded with a software system with all handle functions before shipment. When users need to upgrade the machine functions, they only need to purchase the handle with that function, plug the handle into the machine, and the software system can automatically identify the new handle to realize machine function upgrade.




2D HIFU



3D HIFU



Vmax



Liposonix



Vaginal

Summary of Advantages

🕒 **Procedure Time**
30 - 45 Minutes

❤️ **Full recovery**
Immediate

⊕ **Number of Treatments**
1 every 12 - 18 months

✔️ **Back To Work**
Immediate

🔑 **Software Perloaded**
All functions' software
perloaded already

⚠️ **Risks & Complications**
Possible mild redness,
swelling and tenderness
post treatment

🪡 **Multifunctional**
5 handles optional

📈 **Upgrade**
By inserting new handles only

📄 **Duration Of Results**
12 - 18 months



Before and after



Before

After



Before

After



Before

After



Before

After



Before

After



Before

After